

Anklizer® II Pump Air Walker

NOTE: Instructions cover multiple product families, images may vary.



Figure 1



Figure 2



Figure 3



Figure 4



Figure 5



Figure 6



Anklizer® II Pump
Air Walker Hi Top



Anklizer® II Pump
Air Walker Low Top

NOT MADE WITH
natural rubber latex

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.



WARNING: If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

APPLICATION INSTRUCTIONS

NOTE: High pressure in the bladders provides more support when walking; lower pressure is more comfortable when sitting or reclining. Pressure can be adjusted by inflating/deflating the air bladders or tightening/loosening the device straps.

1. Unfasten straps, open liner, remove pads located inside boot and set aside.
2. While seated, place leg in brace and position heel at the back of the walker. Brace uprights should fall mid-line with the ankle. (Figure 1)
3. Wrap liner around the leg and foot ensuring the ankle pad is positioned over the top of the foot. Secure hook to the liner. (Figure 2)

NOTE: Liner and strap adjustment may be necessary as swelling patterns change throughout the rehabilitation process.

4. Position the uprights parallel to the leg per healthcare provider's recommendation. (Figure 3)
5. Remove plastic sleeves from the uprights to expose hook.
6. Press hook on inside of the upright to the liner. (Figure 4)
7. Beginning with the toe strap, thread straps through links and work up the leg by securing each strap.
8. Alternate directions of the top leg straps and secure to ensure the liner and uprights don't migrate. (Figure 5)

NOTE: Pads may be inserted at either side of the ankle or at heel for additional comfort.

NOTE: Tall ankle walker has three leg straps and the low top ankle walker has one leg strap.

NOTE: The inflation tubes can be re-positioned on the liner or removed if desired.



WARNING: Do not over-inflate the air bladders. Improper inflation may cause skin irritation and/or damage to the device.

9. Inflate air chamber by first turning knob clockwise until knob offers gentle resistance. Press and release air bulb repeatedly until the liner feels stable and supportive on the leg. There should be equal pressure on both sides of the leg, from calf to ankle. (Figure 6)
10. Air can be released from the air chamber by slowly turning the knob counter-clockwise, and halting air release by turning the knob clockwise until the knob offers gentle resistance. To re-inflate, repeat beginning at Step 9.

NOTE: Please refer to sew-in label on device for care instructions.

Please visit our website at www.birdcronin.com for additional product and warranty information or contact a customer care representative at **800-328-1095**.

NOTICE: This device is only one element in the overall treatment program administered by a medical professional. While every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.