

Maternity Sacroiliac Belt

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.

NOTE: Provides comfort when sitting in a car, lounge or office chair. Designed to provide relief from lower back pain and fatigue associated with poor posture, muscle strain, a slipped disc, or back injury.



APPLICATION INSTRUCTIONS

1. Loop the right end of the Maternity Sacroiliac Belt through the metal ring on the 8" or 15" expansion strap. Secure the hook and loop closure.
2. Place the elastic abdominal pad under the natural contour in the center of the abdomen. The wide end should be facing up. The hook and loop closure and sew-in tag will be facing away from the body.
3. Loop and secure the left end of the belt through the remaining metal ring.
4. Adjust tension as needed for comfort and support.

NOTE: Please refer to sew-in label on device for care instructions.

Please visit our website at www.birdcronin.com for additional product information or contact a customer care representative at **800-328-1095**.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

**NOT MADE WITH
natural rubber latex**

