

## B-Cool® Knee Immobilizer

**NOTE:** Instructions cover multiple product families, images may vary.



Figure 1




Figure 2




Figure 3



Figure 4

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

**Intended for single patient use.**



**WARNING:** If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

### APPLICATION INSTRUCTIONS

1. Unfasten the leg straps and unfold the knee immobilizer. (A)
2. Insert one B-Cool® ice pack into each pocket. (Figure 1)
3. Detach and set aside the two 3" hook strips (B) and the two side stays (C) (medial and lateral). (Figure 1)
4. Sit on a flat surface with the affected leg in full extension.
5. Center brace under the leg and align the posterior stays along the back of the leg. (Figure 2)
5. Wrap brace around the leg and hold in place using the two 3" hook strips. (Figure 3)
6. Center a side stay (C) along the inside (medial) length of the leg and press to attach to the brace. Center remaining side stay (C) along the outside (lateral) length and attach. (Figure 3)
7. Wrap and fasten the hook-and-loop straps located directly above and below the knee. Wrap and fasten remaining straps until desired level of fit is reached. (Figure 4)

**NOTE:** Number of straps vary depending on length of immobilizer.

**NOTE:** Universal wrap can be trimmed to reduce excessive amounts of overlap.

**NOTE:** Please refer to sew-in label on device for care instructions.

Please visit our website at [www.birdcronin.com](http://www.birdcronin.com) for additional product information or contact a customer care representative at **800-328-1095**.

**NOTE:** Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

**NOT MADE WITH  
natural rubber latex**