



Bird and Cronin Inc.

Bird & Cronin Inc.

BICRO-LITE™ Ankle Stabilizer



- Designed to help support and protect unstable or injured ankles
- Air and gel bladders are combined to provide compression, support and cold therapy
- Re-engineered shell is lighter in weight but stronger, yet conforms to the ankle with minimal bulk in the shoe
- Anatomically contoured universal shell fits either left or right ankle
- Closed cell foam on the shell provides additional cushioning
- Adjustable heel pad, vertical side straps and swivel straps provide easy application and proper fit
- Regular and small for acute indications
- Trainer/Sport narrow design primarily for chronic indications
- Latex Free

Clamshell Packaging

Reorder No.	Size	Description	Unit
5000 0582	Small	8 1/2" Height	Each
5000 0583	Regular	10 1/2" Height	Each
5000 0584	Trainer/Sport	9" Height	Each