


DRAINAGE BULB HOLDER INSTRUCTIONS FOR USE

Designed to securely hold suction reservoirs during post-operative recovery.



Please read the following instructions carefully and completely before use.
Correct application is required to ensure proper function of the device.

Intended for single patient use.

APPLICATION INSTRUCTIONS

1. Stretch Drainage Bulb Holder around body and press white hook tab to secure.
2. Rotate Drainage Bulb Holder until the locking tabs are lined up with drainage bulbs.
3. Lift blue hook locking tab from holder and thread through loop attached to drainage bulb. Press hook to secure. Repeat as necessary for up to four drainage bulbs.

CARE INSTRUCTIONS

Secure Hook to holder to close. Hand wash with cold water. Air Dry.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.


Please visit our website at BirdCronin.com for additional product information or contact a customer care representative at 800-328-1095.

⚠️ WARNINGS:

If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

DRAINAGE BULB HOLDER INSTRUCTIONS FOR USE

Designed to securely hold suction reservoirs during post-operative recovery.



Please read the following instructions carefully and completely before use.
Correct application is required to ensure proper function of the device.

Intended for single patient use.

APPLICATION INSTRUCTIONS

1. Stretch Drainage Bulb Holder around body and press white hook tab to secure.
2. Rotate Drainage Bulb Holder until the locking tabs are lined up with drainage bulbs.
3. Lift blue hook locking tab from holder and thread through loop attached to drainage bulb. Press hook to secure. Repeat as necessary for up to four drainage bulbs.

CARE INSTRUCTIONS

Secure Hook to holder to close. Hand wash with cold water. Air Dry.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

Please visit our website at BirdCronin.com for additional product information or contact a customer care representative at 800-328-1095.

⚠️ WARNINGS:

If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.