

ABDOMINAL BINDER WITH DRAIN FASTENERS INSTRUCTIONS FOR USE



Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.

APPLICATION INSTRUCTIONS

1. Position the abdominal binder between the lower rib area and hip to cover the abdominal area.
2. Secure the hook and loop closure with enough stretch to provide comfortable support to the abdominal area.
3. Lift blue hook locking tab from binder and thread through loop attached to drainage bulb. Press hook to secure. Repeat as necessary for up to four drainage bulbs.
4. Re-adjust as needed for comfort.

⚠️ WARNINGS:
If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

NOTE: Refer to sew in label for care instructions.

Please visit our website at BirdCronin.com for additional product information or contact a customer care representative at 800-328-1095.

ABDOMINAL BINDER WITH DRAIN FASTENERS INSTRUCTIONS FOR USE



Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.

APPLICATION INSTRUCTIONS

1. Position the abdominal binder between the lower rib area and hip to cover the abdominal area.
2. Secure the hook and loop closure with enough stretch to provide comfortable support to the abdominal area.
3. Lift blue hook locking tab from binder and thread through loop attached to drainage bulb. Press hook to secure. Repeat as necessary for up to four drainage bulbs.
4. Re-adjust as needed for comfort.

⚠️ WARNINGS:
If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

NOTE: Refer to sew in label for care instructions.

Please visit our website at BirdCronin.com for additional product information or contact a customer care representative at 800-328-1095.